

SURFACE DISINFECTANT

READY TO USE

SPORICIDAL DISINFECTION & CLEANING FOR ANY SURFACE

Sychem Control Plus is a patented, proven HLD (High Level Disinfectant). The formulation does not contain any hazardous aldehydes or chlorine generating components.

Sychem Control Plus is highly effective as a cleaner and terminal HLD against bacterial spores, mycobacteria, bacteria, fungi and viruses. The product is designed to eliminate the risk of cross-infection from spores, TB, viruses, fungi and bacteria. Sychem Control Plus does not simply render the microorganisms inert, but has been demonstrated to kill these by disrupting the RNA of the microorganisms. It is safe to use on any surface and is odourless and colourless.

The new version of Sychem Control Plus follows the current Global Harmonisation System (GHS) for labelling and hazard identification, and complies with the REACH 2020/878 (previously 1907/2006) and as amended by EU regulatory 453/2010.

Sychem Control Plus Key Features:

- Sporicidal disinfection and cleaning
- Does not contain chlorine or aldehyde compounds.
- Effective against bacterial spores e.g. C. diff, Bacillus cereus etc
- Effective over a wide range of pathogens including mycobacteria, viruses, bacteria, fungi, yeasts and moulds, including MRSA, Hepatitis B and HIV
- Effective against the difficult non-enveloped type viruses such as Polio, Adenovirus and Norovirus
- Non-fixative
- Non-corrosive - will not damage delicate coatings
- No activation required
- Cost effective
- No fumes
- Tested to EN14885 for disinfection properties



SURFACE DISINFECTANT

Application

READY TO USE

General Surface Cleaning:

For High Level Disinfection, use Control Plus trigger spray.

Spray Control Plus solution onto a cloth/wipe and apply to the surface being treated. Wipe over the whole area.

Allow a minimum contact time of 1 minute before wiping clean or brushing to remove grease, blood and organic debris. Removal of stubborn deposits may be assisted by the use of a stiff bristle brush. Rinse treated surface with clean water if required.

Control Plus solution has demonstrated efficacy in the presence of organic soil contamination and microbiological burden during reuse. The active constituents' concentration of Control Plus solution during its use-life must be verified by the Control Concentration Test Kit or Test Strips.

Sychem Control Plus Microbiological Properties

Active Against	Standards	Contact Time
Spores	EN13704	2 Minutes
	EN14347	1 Minute
Mycobacteria	EN14348	1 Minute
	Mycobacterium terrae	1 Minute
	Mycobacterium tuberculosis	1 Minute
	EN14563	5 Minutes
Bacteria	Mycobacterium avium	5 Minutes
	EN1276, EN13727	1 Minute
	EN13697, EN14561	5 Minutes
Viruses	EN14476	5 Minutes
	Polio, Adeno & Noro	1 Minute
	HIV, Hep B & C	1 Minute
Fungi/Yeasts	EN1275, EN1650	1 Minute
	EN13624, EN14562	5 Minutes

Active ingredients: Bis (3-aminopropyl) dodecylamine
Didecyldimethyl ammonium Chloride.